

# CURIE

## Versatility that Delivers

### COMPACT AND EFFICIENT

Engineered to offer sufficient curing chamber space while maintaining a small footprint.

### SURROUND CURING

UV-lights are strategically positioned to provide a consistent 360° light source for curing and to ensure durability.

### IMPROVED VERSATILITY

365/385/405 nm mix-matrix LEDs for improved rigidity and suitability for for Class 2 bio-materials.



## PRODUCT FEATURES



### **Efficient Curing Space**

Perfect curing chamber space for Dental Clinics and Labs with rigorous workflows.



### **Future-proof**

Equipped with an SD card reader for updating material settings and post-curing times.



### **Accuracy and consistency**

Preloaded optimal curing times database available at the turn of a dial and touch of a button.




### **Premium Performance**

Smart UV-light source management system for maximized energy efficiency and the ideal light strength to prevent any deformation.

# CURIE

- CURIE is equipped with a sophisticated UV-light source management system that enables users to do more in less time and space.
- Produce consistent results thanks to preloaded curing times and easily add new material settings with its built-in SD Card reader.
- CURIE empowers you to deliver premium quality to your stakeholders at affordable prices and in a short time.



# ACKURETTA 3D PRINTING SOLUTION

Dental clinics and labs now have the tools to effortlessly 3D print end-product restorations chairside with increased efficiency, and at a reduced cost. Ackuretta offers a validated all-in-one biocompatible solution that includes hardware, software, finishing solutions, resins, training resources, and documentation that will enable you to complete your prints in a single day.



## 3D Printing Resins

Ackuretta offers a complete and versatile portfolio of CURO resins specially designed for dental applications which deliver high-strength, durable prints for enhanced patient comfort.



## LCD 3D Printers

Same-day chairside dental 3D printers for digital dentistry. Ackuretta's easy-to-use printers support consistent, high-quality prints for dental professionals.

## Proprietary Software

ALPHA AI is the most intuitive dental-focused AI-driven 3D slicing software on the market.

**ALPHA AI**



## Finishing Kit

Enhance your 3D printed dental applications with the Ackuretta Finishing Kit. Clean and cure your prints to ensure safe and ready-to-use outcomes that are professional-grade.

