

# SOL

The perfect **SAME-DAY**  
chairside solution for  
**DIGITAL DENTISTRY.**

**COMPACT, POWERFUL,  
AND ACCURATE,** ready  
to print right out of the box.

Embrace SOL for your  
**DENTAL WORKFLOW**  
today, and start creating  
smiles for the future.



## PRODUCT FEATURES



54 LEDs calibrated for extra curing consistency



150+ Prevalidated Materials




Selective curing Extended accuracy and lifespan



49µm pixel-size Super smooth surface finish



Instant access to resin updates



**PacDent**  
Passion for Excellence

**VOCO**

**Whip Mix**

**DMG**

**NextDent®**

**DENTCA™**  
CAD/CAM DENTURE

**ICTUS**  
DENTAL

**HpDent**

**SHERA**  
IWERKSTOFF  
TECHNOLOGIE

**BEGO**

**ARGEN®**

**NYTE 3D**

**MACKLED**  
RESIN-EXPERTS

**dentona®**  
NEW WAY TO CURE TECHNOLOGY

**SAREMCO**

**keystone**  
industries

# APPLICATIONS



**Dental Model**

3 per print  
54 min



**Temporary Crown**

45 per print  
51 min



**Splint**

3 per print  
2h 07 min



**Denture Base**

3 per print  
1h 51 min




**Denture Teeth**

3 per print  
56 min



**Surgical Guide**

10 per print  
53 min



**IBT**

6 per print  
51 min



**Castable Crown**

45 per print  
1h 25 min



**Impression Tray**

6 per print  
3h 37 min



**Gingiva**

6 per print  
60 min

\* The printing times and pieces per print are obtained using the Large Build Platform and the optimal settings for our CURO resin range.

## ACKURETTA 3D PRINTING SOLUTION

Dental clinics and labs now have the tools to effortlessly 3D print end-product restorations chairside with increased efficiency, and at a reduced cost. Ackuretta offers a validated all-in-one biocompatible solution that includes hardware, software, finishing solutions, resins, training resources, and documentation that will enable you to complete your prints in a single day.



### 3D Printing Resins

Ackuretta offers a complete and versatile portfolio of CURO and FLOW resins specially designed for dental indications which deliver high-strength, durable and biocompatible prints for enhanced patient comfort.




### LCD 3D Printers

Same-day chairside dental 3D printers for digital dentistry. Ackuretta's easy-to-use printers support consistent, high-quality prints for dental professionals.

### Proprietary Software

ALPHA AI is the most intuitive dental-focused AI-driven 3D slicing software on the market.

**ALPHA AI**



### Finishing Kit

Enhance your 3D printed dental applications with the Ackuretta Finishing Kit. Clean and cure your prints to ensure safe and ready-to-use outcomes that are professional-grade.

